


# PULSE CHECK



## AUTOCORRELATOR

Accurate pulsewidth measurement with our family of autocorrelators



"Plug and play", ease of alignment

Femtosecond resolution

Background-free and collinear ACF

Fringe resolved traces

Compact and rugged construction

Pulse duration readout

Computer interface

Laser Diagnostics

Spectral Analysis

Acousto-optics

Non-Linear Optics

Accessories

**Your Partner in Ultrafast**

# PULSE CHECK

## AUTOCORRELATOR

### SPECIAL DELAY LINE

APE PulseCheck has a specially designed, spring loaded linear delay drive which allows you to choose the scan range according to the pulse duration. This offers high accuracy and large duty cycle. Because of the friction-free movement there are scan ranges down to 150 fs. The delay is measured online with high resolution to stabilize scan amplitude and to ensure a linear, calibrated time scale. This also provides the capability to measure interferometric autocorrelation functions.

### ONLINE POSITION MEASUREMENT

The PulseCheck optics is designed using high precision all-reflective optical elements in a very compact way. It comes factory aligned, which means minimal adjustment efforts and easy beam handling for the customer. Minimal dispersion guarantees highest possible resolution. The optics unit incorporates beam shutters for convenient SHG signal check and a servo driven angle tuning of the SHG crystal.


### UNIQUE OPTICAL DESIGN

The PulseCheck offers a simple and continuous adjustment of the interaction angle inside the SHG crystal making it easy to switch between collinear and background-free autocorrelation mode. For lasers with low repetition rate a triggered mode has been implemented to determine the autocorrelation function from the superposition of single pulses.

### MENU SELECTABLE FUNCTIONS AND INTEGRATED DISPLAY

The PulseCheck microprocessor control unit provides menu selectable functions including averaging and data storage. A bright graphical colour display shows the autocorrelation function together with an alphanumerical readout and analog bargraph indication of actual autocorrelation half width. The autocorrelator comes standard with an RS232 serial port.

## AUTOCORRELATOR OPTICS




## DIMENSIONS (in mm)

Control Unit (W\*L\*H):

275\*279\*240 (max.)

Optics Unit:



## SPECIFICATIONS

Version	15	50	150	SM
Scan ranges	150 fs ... 15 ps	500 fs ... 50 ps	1,5 ps ... 150 ps	2,5 ps ... 250 ps <sup>1)</sup>
Delay resolution	< 1 fs	< 2 fs	< 6 fs	10 fs
Min. pulse width	50 fs <sup>2)</sup>	50 fs <sup>2)</sup>	120 fs <sup>2)</sup>	120 fs
Max. pulse width	3,5 ps	12 ps	35 ps	> 60 ps
Scan rate, appr.	13 Hz	10 Hz	7,5 Hz	10 ps/s
	1) Larger scan ranges optional		2) Short-pulse options available	
Linearity	Better 1% of actual scan range			
Sensitivity ( $P_{AV} * P_{PEAK}$ )	PMT: $10^{-4} W^2$ (higher sensitivity optional); PD: $1 W^2$			
Wavelength ranges (others optional)	VIS I		420 ... 550 nm	
	VIS II		540 ... 750 nm	
	NIR		700 ... 1100 nm	
	IR		1000 ... 1600 nm	
	Cross I		360 ... 450 nm + 720 ... 900 nm	
	Cross II		260 ... 320 nm + 730 ... 960 nm	
Input polarization	Horizontal (polarization rotator optional)			
Min. laser repetition rate	Any (optional photodiode detector recommended for repetition rates < 50 kHz)			
Interaction	Non-collinear ... Collinear (intensity and fringe resolved ACF)			
Power supply	95 ... 240 V, 50 ... 60 Hz, 60 W			
Readout	Colour LCD 320 x 240 pixels			
Outputs	Delay: analog		0 ... 10 V	
	Signal: analog		0 ... 10 V	
	RS232 serial interface (parallel interface IEEE488 and USB serial interface optional)			
Input	Trigger: TTL, < 10 kHz			

## OPTIONS

- Spectrometer
- FROG - option
- Additional optics sets
- Fibre input
- Short pulse option
- Triggered SlowScan (recommended for repetition rates < 100 Hz)
- Logarithmic preamplifier
- Enhanced Sensitivity
- Customized wavelength ranges
- Windows control software
- LabView driver
- IEEE488 and USB serial (instead of RS232)
- Input polarization rotator

Distributors  
see APE website